




# OPEN OVAL SHAPE JAWS


With this forceps, you can open a window in the Silicone stent you placed in your patient. This window will allow the ventilation of the safe surrounding bronchi (ex: Right upper Lobe)




The quality of the cut minimize the risk of Granulomas formation

With the same forceps, you may perform different size and shapes of perforation:

- Round,
- Oval,
- Heart shape
- ...



Quality of the cut  
**CLEAN & PRECISE**



Reference:	SCF
<b>Technical Features</b>	
Material:	Stainless Steel
Shaft Length:	180,00 mm
Head Diameter:	8,00 mm
Head Active Length:	10,00 mm

\*too strong material cut may damage the forceps

